


Self test devices – Reading results

COVID-19 Self-Test (Rapid Antigen Test)



Positive result


Two lines appear. One coloured line should be in the control line region (C), and another coloured line should be in the test line region (T).

Two lines, one next to C and one next to T, even faint lines, show the test is positive.



Negative result


One coloured line appears in the control line region (C). No line appears in the test line region (T).



Void result

Control line (C) fails to appear.

Orient Gene




Positive result

The presence of two lines, a control line (C) and a test line (T), even a faint line, within the result window indicates the test is positive.



Negative result


The presence of only a control line (C) with the result window indicates the test is negative.



Void result

If the control line (C) is not visible within the result window after performing the test, the result is considered void. If you get void result, you must perform a new test.

ACON Flowflex




Positive result

Two lines appear. One coloured line should be in the control line region (C), and another coloured line should be in the test line region (T). Two lines, one next to C and one next to T, even faint lines, shows the test is positive.



Negative result

One coloured line appears in the control line region (C). No line appears in the test line region (T).



Void result

Control line (C) fails to appear.